Accepted Manuscript

CHROMATOGRAPHY OF OXYSTEROLS

Irundika H.K. Dias, Steven R. Wilson, Hanne Roberg-Larsen

PII: S0300-9084(18)30127-5

DOI: 10.1016/j.biochi.2018.05.004

Reference: BIOCHI 5416

To appear in: *Biochimie*

Received Date: 26 February 2018

Accepted Date: 4 May 2018

Please cite this article as: I.H.K. Dias, S.R. Wilson, H. Roberg-Larsen, CHROMATOGRAPHY OF OXYSTEROLS, *Biochimie* (2018), doi: 10.1016/j.biochi.2018.05.004.

This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.



Abstract

Oxysterols play important roles in development and diseases, but can be highly challenging to analyze. To ensure satisfactory measurements, oxysterols must typically be separated with chromatography prior to detection. Here, we will devote attention to the chromatography of oxysterols, focusing on gas chromatography and liquid chromatography. We will present the role of stationary phases, mobile phases, and dimensions and geometries of particles/columns. We discuss how these parameters may affect the chromatography, regarding factors such as speed and resolution. Finally, we present some less explored avenues for separation of oxysterols.

- 1 CHROMATOGRAPHY OF OXYSTEROLS
- 2 Irundika HK Dias^a, Steven R Wilson^b and Hanne Roberg-Larsen^{b*}.
- ^aAston Medical Research Institute, Aston Medical School, Aston University, Birmingham, B4
- 4 7ET, UK
- 5 ^bDepartment of Chemistry, University of Oslo, PO box 1033, Blindern, NO-0315 Oslo,
- 6 Norway
- 7
- 8 * Corresponding author: Hanne Roberg-Larsen, Department of Chemistry, University of Oslo.
- 9 PO box 1033 Blindern, NO-0315 Oslo, Norway. Email: <u>hanne.roberg-larsen@kjemi.uio.no</u>.
- 10 Phone +47 22 85 55 87
- 11



24 Oxysterols; Chromatography; Separation; Lipids; Sterols

25 Abstract

26 Oxysterols play important roles in development and diseases, but can be highly challenging to 27 analyze. To ensure satisfactory measurements, oxysterols must typically be separated with 28 chromatography prior to detection. Here, we will devote attention to the chromatography of 29 oxysterols, focusing on gas chromatography and liquid chromatography. We will present the 30 role of stationary phases, mobile phases, and dimensions and geometries of particles/columns. 31 We discuss how these parameters may affect the chromatography, regarding factors such as 32 speed and resolution. Finally, we present some less explored avenues for separation of 33 oxysterols.

34	
35	
36	
37	
38	
39	
40	
41	
42	
43	
44	
45	
46	
47	
48	
49	
50	
51	
52	
53	
54	
55	
56	
57	

58

59 **1. Introduction**

Oxysterols are a group of lipids that receive considerable attention due to the unraveling of 60 61 their roles in numerous diseases and development [1-5], and are established biomarkers for 62 e.g. Niemann-Pick disease (NPD) [6]. Quality measuring tools must be employed to 63 understand the roles oxysterols play in development, diseases and conditions. However, the 64 measurement of oxysterols can be highly challenging. Some reasons are that oxysterols may be present at low concentrations, in limited samples. In general, such issues can often be 65 66 solved by using highly sensitive mass spectrometry (MS) techniques. However, many 67 oxysterols are not "ideal" for MS analysis, as they can be difficult to ionize; ionization is key 68 requisite when using electrospray ionization (ESI), a most common interface of MS. But 69 perhaps equally important, oxysterols are often highly similar compounds, e.g. present as 70 isomers with similar MS fragmentation profiles, making selective determinations a significant 71 challenge. Thus, oxysterols require particular care regarding pre-MS steps. A key step is to 72 ensure high quality chromatographic separations, for e.g. resolving isomers and achieving 73 precise measurements.

74 In this review, we will focus on the chromatography of oxysterols. In particular, we will 75 discuss separations of oxysterols using gas chromatography (GC) and liquid chromatography 76 (LC), giving attention to the speed, resolution and sensitivity of oxysterol separations using 77 these techniques. Although we acknowledge the great improvements made in separation 78 instrumentation over the years (rapid injection systems, low void volume connections, MS 79 resolution etc.), we focus here on fundamental separation conditions, e.g. column materials, 80 stationary phases, mobile phases (MP) and particle geometries. The chromatography of native 81 and derivatized oxysterols ("charge-tagged" for improved MS sensitivity) will be discussed. 82 Finally, we will take a look at some less employed approaches for separations, which may 83 have future roles in oxysterol separations.

84

85 2. Gas Chromatography and oxysterols: "Never change a winning team"?

GC is a technique in which compounds are separated in meter-scale columns with inner diameters well below 1 mm. Compounds are separated by having unequal retention factors (time spent on the column walls/ time spent in a gaseous MP). The stationary phase is typically a polymer coating around 0.25 μ m in thickness. GC can provide excellent resolution and is simply coupled with MS, typically via election ionization (EI) interfaces. Also, GC-EI-MS does not suffer from suppression effects to the same degree as ESI (the common MS

92 interface with LC) [8, 9]. Reduced suppression from other compounds lessens the need for 93 analyte-specific internal standards. GC has been a workhorse for analysis of sterols for well 94 over 50 years [10-12]. For the last couple of decades, a method described by Ulf Diczfalusy 95 and co-workers for oxysterol analysis has been highly influential [13]. For measurements of 96 the analytes in human plasma, the authors separated 7α - and 7β -hydroxycholesterol, 7-97 oxocholesterol, cholesterol- 5α , 6α -epoxide, cholesterol- 5β , 6β -epoxide, cholestane- 3β , 5α , 6β -98 triol, 24-hydroxycholesterol, 25-hydroxycholesterol, and 27-hydroxycholesterol and used MS 99 for detection (Figure 1). The analytes were derivatized with trimethylsilyl (TMS) to enhance 100 volatility as required by GC. The authors employed a 30 meter long column, featuring a HP-101 5MSM ((5%-phenyl)-methylpolysiloxane) film, which is described as a non polar "general 102 purpose" stationary phase, where separations are mainly based on differences in boiling 103 points. A temperature gradient was employed and the MP (helium) was held at 0.8 mL/min 104 (35 cm/second), which is close to the optimal column efficiency for the MP and column inner 105 diameter (ID) (0.25 mm). Using these conditions, the authors obtained peaks with widths of 106 about 15 seconds, with all the analytes being detected within 19 minutes.

107 This method has been cited several hundred times, and more importanty, has been 108 reproduced/re-used in a substantial number of studies. Quite remarkably, few major 109 significant modifications to this GC method have been reported. Some researchers have 110 reduced the analysis time; for determining serum cholestane- 3β , 5α , 6β -triol as a biomarker for 111 Niemann-Pick type C disease (NPC) diagnosis, Kannenberg et al. performed analyses below 112 10 minutes per sample. The authors employed a trifluoropropylmethyl polysiloxane phase, a 113 mid-polar stationary phase, which is promoted as being ideal for separating positional isomers 114 [14]. A most notable reduction of analysis time was shown by Maria T. Rodriguez-Estrada 115 and co-workers, who performed analysis of cholesterol oxidation products using "fast gas 116 chromatography" [15]; a shorter (10 m) column with a narrower ID (0.1 mm) and thinner 117 stationary phase film (0.1 μ m thickness) allowed for separation (resolution (Rs)>1.2) of 7 α hydroxycholesterol, 19-hydroxycholesterol, 7β -hydroxycholesterol, β -epoxycholesterol, α -118 119 epoxycholesterol, cholestanetriol, 25-hydroxycholesterol, 7-KC in 3.5 minutes (Figure 2). It 120 should be noted that a narrow ID and thin film allow for improved efficiency, which can be 121 essential for obtaining resolution when the column length is decreased (reduces analysis time, 122 but affects the separation). The fast analysis time of this method is roughly 5-10 times faster 123 than that described in most related papers. The flow rate was 0.41 mL/minute, which 124 corresponds to a linear velocity of 43 cm/second (narrow columns can be operated at higher 125 linear velocities without dramatic decreases in efficiency). However, substantial efforts

126 towards gas chromatographic improvement have in general not been prioritized in oxysterol 127 analysis. This is perhaps quite reasonable, as GC has to a large degree reached a level of 128 maturity allowing for stable analysis, especially when analysis time is not a major concern. 129 Also, GC-MS instruments are typically less expensive than LC-MS instruments.

130

131 **3. Liquid chromatography and oxysterols: a broad range of approaches**

132 Briefly, LC is a technique in which compounds are separated in a centimeter-scale column 133 with an inner diameter typically between 1-4 mm. The MPs are liquids, and the stationary 134 phases are often attached to particles with sizes typically between 1.7 and 5 µm diameters. 135 Although GC can provide greater plate numbers per column, LC provides greater plate 136 numbers per column length. Moreover, derivatization is not an absolute requirement for all 137 applications. Regarding oxysterols, LC-ESI-MS is generally more sensitive than GC-MS 138 variants [16]. Oxysterol separations by LC have been undertaken using either normal-phase 139 (NP) LC or reversed-phase (RP) LC.

140

141 *3.1 NPLC: a lost cause?*

142 In NPLC, molecules are adsorbed to the surfaces of silica particles (unmodified or featuring polar chemically bonded phases) and are eluted with non-polar solvents such as hexane or 143 144 heptane. Analytes engage with the stationary phase via hydrogen bonding and dipole 145 interactions. There are some examples of determination of oxysterols using NPLC coupled to 146 UV detection with good chromatographic resolution between the isomers [17, 18]. However, 147 NP MP solvents show low conductivity, surface tension, and lack ability to donate or accept a 148 proton to give analytes charge (and hence sensitivity) for ESI-MS detection (a "gold-149 standard" in the analysis of fluids). On the other hand, this can be partially overcome by 150 addition of polar solvents such as 2-propanol [13] or methanol [14]. But in addition, NP is 151 generally perceived as having lower reproducibility and predictability than RPLC. Normal 152 phase and reverse phase separations have been compared regarding oxysterols. Careri et al. 153 compared chromatographic separations of cholesterol and five oxysterols using NPLC 154 (Nucleosil 5-CN column) and RPLC (C18, Ultracarb ODS (20) column). RPLC (associated 155 with more robustness) provided superior sensitivity for all the compounds measured [19]. 156 RPLC is the key mode when analyzing oxysterols in complex biological samples with MS 157 detection. However, we will discuss promising separation techniques related to NPLC in the 158 final section.

159

160 3.2 RPLC: a workhorse with possibilities for different selectivity

161 RP stationary phases are hydrophobic and are typically alkyl chains. The oxysterols, both 162 native and derivatized are most commonly separated using octadecyl alkyl chain (C18)-163 bonded silica stationary phases, followed by C8-bonded silica columns (Tables 1 and 2). The 164 analytes are expected to elute according to hydrophobicity (although RPLC is far more 165 complex than commonly perceived [20]). With RPLC, the more polar side chain oxygenated 166 oxysterols elute before ring-oxidized sterols followed by more non-polar sterols [19, 21]). 167 C18 columns generally provide the same retention order for oxysterols (when using similar 168 MPs) regardless of the derivatization reagent (derivatization in LC is discussed in some more 169 detail below). This implies that selectivity differences are typically not attributed to the 170 derivatization reagent. However, RPLC can have somewhat surprisingly degrees of selectivity 171 Shan et al. compared oxysterol separations regarding two solvent systems, options. 172 acetonitrile:water and methanol:water, on C8 and C18 columns [22]. They showed that even 173 with same MP the chromatographic mobility and selectivity between C8 and C18 columns are 174 considerably different. The C8 column was able to resolve several oxysterol pairs, including 175 7α - and 7β -hydroxycholesterol, which were inseparable on the C18 column using a similar 176 MP. Compared to using methanol:water, using acetonitrile:water with a C18 column provided 177 improved resolution of oxysterol pairs 27-hydroxycholesterol/ 3β , 5α , 6β -triol, 24R-178 hydroxycholesterol/20α-hydroxycholesterol, 7α -hydroxycholesterol/ 7β -hydroxycholesterol, and 7-ketocholesterol/3β-OH-6-one. Roberg-Larsen et al. observed that employing an 179 180 acetonitrile-based MP and C18 allowed separation of 20-hydroxycholesterol and 27-181 hydroxycholesterol, while separation of 27-hydroxycholesterol and 24S-hydroxycholesterol 182 was not possible with the same conditions [23]. The opposite was observed when using 183 methanol-based MP (in the final section we will discuss an approach that may allow these two 184 systems to be combined). Roberg-Larsen et al. has observed highly similar RPLC oxysterol 185 separations in microbore LC, capillary LC and nano LC [23-25] (selectivity is rarely affected by column diameter). 186

187

188 24S-hydroxycholesterol and 25-hydroxycholesterol can also be challenging to separate. 189 Debarber et al. separated 24R- and 24S-hydroxycholesterol by modifying a method by 190 Burkard and coworkers [26] using a methanol:acetonitrile:water MP (45:40:35) and a column 191 temperature of 55 $^{\circ}$ C [27]. Changing to methanol/acetonitrile/water (14:0.6:1) and a column 192 temperature of 10 $^{\circ}$ C same authors demonstrated separation of 24-hydroxycholesterol from 193 25-hydroxycholesterol within short (6.5 min) time. The later method reversed the retention

order of two oxysterols to 24-hydroxycholesterol followed by 25-hydroxycholesterol.
However, the later method did not baseline separate the two oxysterols.

196

197 Great efforts have been made to ensure that the analytes do not interact with the particles in 198 which the stationary phase is attached, as e.g. silica particles cause secondary interactions and 199 may perturb the separation (e.g. cause band broadening and tailing). Avoiding interactions 200 with the particles can be done by efficient rest-silanization and adding functional groups at the 201 trunk of the main stationary phase for steric hindrance. In addition, the carbon loading of the 202 particles is crucial for oxysterol separations. For instance, otherwise high quality columns 203 which featured lower carbon loads were unable to provide selective separation of side-chained 204 oxysterols [25].

205

To improve chromatographic efficiency further, the size of the particles may be reduced; 206 207 today sub-µm particles are common, while 3-5 µm particles were standard about a decade ago. The Hypersil GOLDTM column is a familiar column in oxysterol RPLC, both in UHPLC 208 209 (sub 2 µm particles) and regular HPLC (3 µm particles), typically with a 2.1 mm ID. The Hypersil GOLDTM columns are endcapped silica-based columns, with a high hydrophobicity 210 211 and medium shape selectivity and polar surface activity. (For column classification see [28] 212 and [29]) This column material has been used for chromatographing both Girard P- and 213 Picolinyl ester derivatized (PED) oxysterols [30-39]. This column seems to not be compatible 214 with Girard T derivates (Rs >1.1, data not published). However, sub-µm particles cause 215 higher back-pressures. Therefore, a highly attractive alternative has been the use of core shell 216 particles. Core shell particles have a solid core and porous shell, that gives high efficiency and 217 fast separations with low back pressure compared to traditionally porous particles [40]. Core 218 shell particles provide similar efficiencies to sub-um particles. A well-known example of core 219 shell particles for oxysterol analysis is that by McDonald et al. who chromatographed 62 220 different sterols, oxysterols and secosteroids from human plasma using two different LC and 221 one GC method. Both the LC methods used core shell particles, with the side-chain 222 oxysterols eluting in 7.5min (total run time 12min) [41].

A notable exception from using octyl chain stationary phases is by Silke Matysik and coworkers who employed a biphenyl phase [42]. Biphenyl stationary phases typically provide increased retention and can have a different selectivity compared to traditional C18/RP phases, as it can provide both π – π interactions and higher hydrogen bonding capacities [43]. In addition to featuring a different phase, the column employed was packed with core shell

particles. The work of Silke Matysik and co-workers demonstrate a quick separation of N,N dimethylglycine (DMG) derivatized oxysterols (8 oxysterols in 8 minutes, see Figure 3)

230

231 3.3 Effect of derivatization on oxysterol chromatography

232 Derivatization of oxysterols is used to enhance sensitivity in MS detection by incorporation of 233 a charge group into the oxysterol. Derivatization of oxysterols can be used with both ESI-MS 234 and atmospheric pressure chemical ionization (APCI)-MS. In addition to enabling enhanced 235 sensitivity, derivatization can also make the sterols more soluble in MPs commonly used in 236 RPLC. Derivatization can also prevent adsorption of the hydrophobic sterols on narrow ID 237 fused silica tubing in use in sensitive nano LC-based systems [7]. Other benefits of 238 derivatization are more easily interpretable MS2 spectra, as fragmentation of the derivatized 239 group usually gives more specific fragmentation [7].

240

A variety of different derivatization reactions for sterols exist, and most common used ones 241 242 for oxysterols are Girard P and T reagent, picolinyl acid and DMG (for end product structures 243 see Figure 4 and a recent review by Yuqin Wang and William J. Griffiths' group summarizes 244 the details for all the most common derivatization reactions [4]. While DMG is mostly used in 245 the context of NPC disease [44-46], Girard P and T is used in the context of neurologic [31] 246 or metabolomic [33] diseases and cancer [23, 24]. All these derivatization reactions are 247 targeting the hydroxyl group. An alternative is to use click-chemistry, to target the double 248 bond between the C5 and the C6 in the sterol structure, e.g. by thiol-ene click-chemistry 249 tagging using a photoinitiator [47]. The click-chemistry generates heteroatom links (C-X-C) 250 and reaction rates can be quick (< 1 minutes) when using a microflow reaction cell.

251

252 Regarding chromatographic performance, there are small differences in the behavior of the derivatized or native oxysterols. Cha et al. [48] has analyzed both native oxysterols from 253 serum samples as silver adducts and picolinyl ester derivatized (PED) oxysterols from CFS. 254 255 Although the chromatograms look very different regarding analysis time, the analysis is 256 performed on two different reversed phase columns; An ACE C18 (3 µm, 150 mm x 2.1 mm 257 ID) and a Kinetex C18 (2.6 µm, 100 mm x 2.1 mm ID, core shell). Although these columns 258 has approximately the same hydrophobicity, they have different shape selectivities and polar 259 surface activities [28] and most importantly, different solid supports (fully porous vs. core-260 shell). It would be interesting to compare the separation of the PED with McDonald et al. 261 [41], and the native oxysterol separation with Roberg-Larsen et al [24], which both uses the

same columns on native and Girard T derivatized oxysterol, respectively. McDonald's
oxysterol-ammonium adducts shows similar chromatography and elutes in the same retention
window as Cha's PED, while Roberg-Larsen's Girard T derivates elutes in the same retention
window as Cha's native sterol.

266

267 3.4 Dimensions and sensitivity

The most popular column dimension used in oxysterol analysis is the 2.1 mm ID format. 268 269 Detection limits for native oxysterols and all the types of derivates are in low ng/mL, suitable 270 for analyzing oxysterols in plasma. However, most of the applications use more than 50 µL 271 plasma or serum in their sample preparation. The relatively high sample volume for the other methods can be challenging if the sample sizes are small, e.g. plasma from mouse and rats. 272 273 Exception is the method from Honda et al. [36] and Xu et al. [38], which both used only 5 µL and picolinyl ester derivatization. Sensitivity in the same range has been achieved with 274 275 Girard T derivatization (in cell sample) using narrow bore columns [23, 24]. In general, the 276 sensitivity will depend on both the efficiency of the sample preparation and column 277 dimensions. The 2.1 mm ID columns with small particles (e.g. > 2 μ m), provide high 278 efficiency separations, but more narrow columns (µm-scale IDs, e.g. nano LC and capillary 279 LC) can be employed when the goal is to enhance sensitivity [23, 24, 49].

- 280
- 281

282 4. Unknown Pleasures? Alternative separation approaches for oxysterols

In addition to conventional LC and GC, there are a number of other separation approachesthat are less explored regarding oxysterols.

285 2D GC [50] means to couple two different GC columns in a single system, to enhance 286 chromatographic resolution. The two columns must have different selectivity, and are 287 connected via a modulator. Fractions elute from the first (usually long) column, and are subsequently chromatographed on a second (usually short) column. A large number of 288 289 chromatograms are generated during an analysis, and dedicated software assembles these into 290 a 2D plot (resemblance of a 2D gel). The combined resolution is in theory the product of the 291 peak capacity of the two columns (in practice, this number is lower). This approach is used in 292 e.g. food and gasoline analysis, but has also been used for mapping sterols [51] (Figure 5). 293 However, the approach is not commonplace, but is commercially available from a number of 294 manufacturers.

295

296 2D LC is a similar variant to 2D GC, where two separation columns are coupled, for example 297 hydrophilic interaction chromatography (HILIC) and reverse phase (RP) LC. 2D LC has been 298 used for lipid analysis ([52]), but not for sterols (to the authors' knowledge). It is worth 299 mentioning here that HILIC highly related to NPLC, has an acceptable stability and is highly 300 MS compatible [53]. Hence, it could be interesting to see if this phase would have promise for 301 oxysterol separations. Although 2D LC can provide very high resolutions, it is arguably less 302 straightforward to operate than 2D GC, as different LC columns are often not compatible with 303 each other's preferred MP solvents. However, it could be interesting so see if 2D LC could 304 fully resolve side chain-hydroxylated oxysterol isomers, by combining methanol:water and 305 acetonitrile:water LC separations in a joint system.

306 Capillary electrophoresis (CE) and related techniques are characterized by an electric field 307 applied across an open tube/column in which the separation takes place. Compounds are 308 separated by charge and hydrodynamic radius, often with unprecedented resolutions. CE has 309 been used for sterol analysis, using organic solvents (non-aqueous CE = NACE) [54]. Since 310 the approach does not require a solvent pump, it is highly suited for miniaturization/chip 311 separations. However, it remains to be a more technically challenging technique compared to 312 LC and GC.

313 Open tubular columns (not filled with a particles) are typically used for GC and CE (and 314 related techniques), but are rarely used in LC. However open tubular LC (OTLC) can provide 315 for excellent chromatography and sensitivity. Such columns are typically 10 µm ID, featuring 316 a stationary phase attached to the inner walls, as in GC. OTLC has been demonstrated 317 regarding oxysterols, and Vehus et al. [49] achieved detection limits of 25 attograms (Girard 318 T derivatized 25-hydroxycholesterol). For comparison, previous high sensitivity methods 319 have achieved detection limits in the femtogram range [23, 36]. OTLC is predicted to have a 320 significant role in tomorrow's liquid separations [55]. However, as with the other techniques 321 presented in this section, it has larger technical challenges, where routine labs cannot be expected to have patience for. This may be resolved when commercial OTLC 322 323 products/systems become available, although these will perhaps be primarily used for 324 applications with very limited amounts of sample. In addition to the techniques described 325 here, there are other approaches that are rather unexplored regarding oxysterol analysis. For 326 instance, supercritical fluid chromatography (SFC) may be an interesting and useful approach, 327 as SFC is associated with speed and ability to separate isomers. SFC has previously been 328 demonstrated with other sterols and related compounds [56, 57].

329

5. Conclusions

331 Oxysterols can be challenging to separate, and some oxysterol pairs such as 24 R/S -332 hydroxycholesterol and 25-hydroxycholesterol, 7α - and 7β -OHC need particularly careful 333 attention to chromatographic separation. Indeed, mass spectrometry can offer an additional 334 level of resolution by differentiating co-eluting compounds by mass and selecting specific ion 335 pairs, e.g. with multiple reaction monitoring methods (MRM). However, quality oxysterol 336 analysis needs quality separations. Regarding oxysterols, LC is becoming increasingly used 337 and developed compared to GC. Newer types of solid support (e.g. core-shell) and stationary 338 phases (e.g. biphenyl) should be further explored for more time efficient separation. 339 Sensitivity is good enough for native oxysterols in serum/plasma if sample sizes are ample 340 $(>100 \ \mu L)$, but the inner diameter of the column can be modified to obtain sensitivity gains. 341 Since chromatography is an important aspect in the analysis of oxysterols, we encourage 342 readers to provide details on their chromatographic methods and challenges, to set the stage 343 for faster and more efficient analyses in the future.

- 344
- 345

346 6. Acknowledgements

347

348 HKID acknowledges support from Kidney Research Foundation PDF3/2014 and Alzheimer's
349 Research UK network grants. This work has been performed within DIATECH@UIO,
350 strategic research initative at the Faculty of Mathematics and Natural Science, University of
351 Oslo, Norway (HRL and SRW).

- 352
- 353
- 354

REF	Derivatization	Column	Dimensions	Particle	Pore	Surface	Carbon load	Mobile phases and	Detection	Run	Analytes in retention order
			L x ID (mm)	size	size	area	(%)	Temperature	limits	time	
				(µm)	(Å)					(min)	
[58]	N-4-(N,N-	Acquity UPLC CSH TM C ₁₈	100 x 2.1	1.7	130	185	15	Formic acid (FA)	Low nM	13	22R-OHC, 27-OHC,
	dimethylamino)phenyl							in			25-OHC, 24S-OHC,
	carbamates							H ₂ O/MeOH/ACN,			7β-ОНС, 5β,6β-
						Ċ		70 °C			epoxycholestanol, 5a,6a-
											epoxycholestanol,
											desmosterol, 7-
											dehydrocholesterol,
						1					lathosterol, cholesterol,
					\mathcal{N}						cholestenol
[44]	N,N-dimethylglycine	BetaSil C ₁₈	100x 2.1	5	100		20	Trichloroacetic	2 ng/mL	10	3β,5α,6β triol, 7-keto-OHC
								acid (TCA)/acetic			
								acid (AA) in			
								H ₂ O/ACN			
[46]	N,N-dimethylglycine	Betasil C ₁₈	100x 2.1	5	100		20	TCA/AA in	2 ng/mL		3β,5α,6β triol, 7-keto-OHC
								H ₂ O/ACN			
[42]	N,N-dimethylglycine	Kintex [™] Biphenyl	50 x 2.1	2.6	100	200	11	FA/Ammonium	1 ng/mL	8	25-OHC, 24S-OHC, 27-
				Core shell				acetate in			OHC, 4β -OHC, 7α -OHC
								H ₂ O/MeOH/ACN,			7β-OHC, 7-keto-OHC,
								30 °C			$3\beta,5\alpha,6\beta$ triol
[59]	N,N-dimethylglycine	Gemini-NX TM C ₁₈	100x2	3	110	375	14	Ammonium	0.08-0.8	15	3β,5α,6β triol, 7-keto-OHC
								formate in	ng/mL		
		Y						H ₂ O/ACN			
[35-37, 60]	Picolinyl ester	Hypersil GOLD C ₁₈	150 x 2.1	3	175	220	10	FA in	2-10 fg on	40	24S-25-epoxy-OHC, 22R-
								H ₂ O/MeOH/ACN,	column		ОНС, 24S-ОНС, 25-ОНС
								40°C			27-ΟΗC, 7α-ΟΗC, 4β-

Table 1 Chromatographic conditions for derivatized oxysterol separation

											OHC,
[61]	Picolinyl ester	Acquity UPLC BEH C ₁₈	100 x 2.1	1.7	130	185	17.7	FA in H ₂ O/MeOH/ACN, 35 °C	2 ng/mL	11	4α-ΟΗC, 4β-ΟΗC
[38]	Picolinyl ester	Hypersil GOLD C ₁₈	50 x 2.1	1.9	175	220		FA in H ₂ O/ACN , 25°C	5 ng/mL	16	24S-OHC/25-OHC*, 27- OHC/7α-OHC/7βOHC*, 4α-OHC, 4β-OHC, cholesterol
[39]	Picolinyl ester	Hypersil GOLD C ₁₈	50 x 2.1	1.9	175	220	10	AA in H ₂ O/ACN, 40 °C	5 ng/mL	15	4α-ΟΗC, 4β-ΟΗC
[48]	Picolinly ester	Kinetex C ₁₈	100x 2.1	2.6	100	200	12	FA in H ₂ O/MeOH, 25 °C	0.5-5 ng/mL	10	24S-OHC, 25-OHC, 27- OHC
[30-34]	Girard P	Hypersil GOLD C ₁₈	50 x 2.1	1.9	175	220	10	FA in H ₂ O/MeOH/ACN		17	24S-OHC, 25-OHC, 27- OHC, 7β-OHC, 7-O- OHC,7α-OHC, 6-OHC
[62]	Girard P	Kinetex C ₁₈	50 x 2.1	1.7	100	200	12	FA in H ₂ O/MeOH/ACN		17	24S-OHC, 25-OHC, 27- OHC, 7β-OHC, 7-O-OHC, 7α-OHC, 6-OHC
[25]	Girard T	ACE C ₁₈	150 x 1	3	300	100	9	FA in H ₂ O/ACN, 40 °C	0.2 nM	20	25-OHC, 24S-OHC, 20α- OHC, 22S-OHC
[23]	Girard T	ACE C ₁₈	150 x 0.1	3	300	100	9	FA in H ₂ O/MeOH	23 pM	40	22R-OHC, 24S-OHC, 25- OHC, 27-OHC, 22S-OHC
[24]	Girard T	ACE C ₁₈	150 x 0.3	3	300	100	9	FA in H ₂ O/MeOH	25 pM	35	22R-OHC, 24S-OHC,25- OHC, 27-OHC, 22S-OHC
[45]	Dimethylaminobutyrate ester	Phenomenex Synergi fusion C ₁₈	50 x 2.1	4	100	475	12	FA + ammomium formate in H ₂ O/ACN	0.5 ng/mL	6	3β,5α,6β triol, 7-keto-OHC

*Coelution

Table 2: Chromatographic conditions for native oxysterol separation

								X		
REF	Column	Dimensions	Particle	Pore	Surface	Carbon	Mobile phases and	Detection	Run time	Analytes in retention order
		L x ID (mm)	size	size	area	load (%)	Temperature	limits	(min)	
			(µm)	(Å)	(m2/g)		0-			
[63]	Zorbax Eclipse Plus C18	150 x 2.1	3.5	95	160	9	Ammonium acetate in H ₂ O/MeOH/ 30°C		30	24-OHC, 27-OHC, desmosterol, cholesterol, lanosterol, cholestanol stigmasterol campesterol, β- sitosterol, sitostanol
[26]	Nucleosil C18 HD	125 x 2	5	120	200	11	Ammonium acetate in MeOH/ACN/ H ₂ O	25 ng/ml	35	24S-OHC, 27-OHC
[19]	Nucleosil 5-CN	250 x 2	5	100	350	5	Heptane/Propan-2-ol	16 ng	20	Cholesterol, $5,6\alpha$ -EP, 25- OHC, 7-keto, 7 β -OHC and $3\beta,5\alpha,6\beta$ triol
	Ultracarb ODS (20) C18	250 x 2	5	90	370	22	MeOH/ACN	4 ng	20	25-OHC, 3β,5α,6β triol , 7β-OHC, 7-keto, 5,6α- epoxy-OHC, cholesterol
[27]	BetaBasic C18	250 x 2.1	5	150	200	13	Ammonium acetate in MeOH/ACN/H ₂ O 10°C	30 ng	30	25-OHC, 24-OHC
[64]	Synergi Hydro	250 x 2	4	80	475	19	MeOH/ACN/H ₂ O 30°C	0.1-0.4 ng/ml	25	3β,5α,6β triol , 7α-OHC, 7β-OHC, 7-keto, β-epoxy- OHC, α-epoxy, 6-keto
[65]	Supleco Ascentis®MS (C8)	100 × 2.1	3	100	450	15	ACN/H ₂ O/Ammonium acetate	4 ng/ml	7	4β-ОНС
[66]	Chromolith SpeedRod RP- 18e monolithic	50 × 4.6	2		250	18	MeOH/H ₂ O	0.1 ng/ml	7	cholestane 3β,5α,6β-triol, 7- α/β-hydroxycholesterol, 5,6- β-epox-OHC, 5,6-α-epoxy- OHC, 7-ketocholesterol,

										cholesterol.
[67]	Nucleosil HD, C18	250 × 4.6	3	100	350	20	FA in MeOH/ H ₂ O	10 ng/ml	14	270НС
[68]	LiChrosorb RP-18	250 ×4	5	100	300	17	MeOH/ACN 30°C	0.2 ng	16	25-OH, 3β,5α,6β triol, 7β- OH, 7-keto,5,6α –epoxy- OHC, cholesterol
[41]	Kinetex C 18	150×2.1	2.6	100	200	12	Ammonium acetate in ACN//IPA	1 ng/ml	12	60 analytes
[69]	ACQUITY UPLCTMBEH C18	150 × 2.1	1.7	130	185	18	FA in MeOH/H ₂ O 40°C	54 pg/ml	20	24-OHC, 25-OHC 7-OHC, 4β-OHC and 7-keto cholesterol
[21]	Shimadzu Shim-pack ODS	100 × 3	2.2	8nm	470	20	H ₂ O/ ACN 50°C		16	24(S)-OHC, 25-OHC, 27- OHC, 7α, 7β, 4α-,5,6β- epoxy-OHC, 5,6α-epoxy- OHC, 4β-OHC, cholesterol
[70]	ODS AQ C18	150 × 4	5	120	330	14	MeOH/ACN/H2O	100 ng/ml	30	25-OHC, cholestane-3β-5α- 6β-triol, 7β-OHC, 7- ketocholesterol, 5,6α-epoxy- OHC, cholesterol.
[71]	Aquasil C18	250 x 4.6	5	100	310	12	ACN/MeOH 25°C	0.5ng	19	7α-, 7β-, 25-OHC, 7-keto, 3β,5α,6β triol, α-epoxy, β- epoxy
[72]	Nova Pack CN HP	300 x 3.9	4	60	120	3	n-Hexane-2-Propanol 32°C	6-70 ng/ml	30	19-OHC, cholesterol, 20 α- OHC, 22(R)-OHC, 24(S) – OHC, 22(S)-OHC, 25-OHC, 5,6 α-epoxy-OHC, 5,6 β- epoxy-OHC, 25(R)-OHC, 7-ketocholesterol, 7β-OHC, 7α-OHC
[73]	NUCLEOSIL® C18	100 x 4	5	100	350	15	FA in MeOH/H ₂ O/2- propanol	5-135 pg/ml	45	24-OHC, 25-OHC, 27- OHC, 7β-OHC, 7- ketocholesterol
[5]	Supelcosil LC-18-S	250 x 4.6	5	120	170	11	FA in MeOH/H ₂ O	3.2 ng/ml	45	21 analytes
			R)					

FIGURE LABELS

Figure 1. GC-MS performance of the method by Diczfalusy and co-workers [13]. Broken lines are unlabeled compounds, and solid lines are deuterated internal standards. Compounds separated (plasma sample) are: I. 7α -hydroxycholesterol, II. 7β -hydroxycholesterol, III. cholesterol- 5α , 6α -epoxide, IV. cholesterol- 5β , 6β -epoxide, V. cholestane- 3β , 5α , 6β -triol, VI. 24-hydroxycholesterol, VII. 25-hydroxycholesterol, VIII. 7-oxocholesterol, IX. 27hydroxycholesterol. All compounds were derivatized with TMS. Reprinted with permission.

Figure 2. GC-MS performance of the method by M T Rodriguez-Estrada and co-workers [15]. The total ion current chromatogram shows a fast GC-MS separation of 1. 7α -hydroxycholesterol, 2. 19-hydroxycholesterol, 3. 7β -hydroxycholesterol, 4. β -epoxycholesterol, 5. α -epoxycholesterol, 6, cholestanetriol, 7. 25-hydroxycholesterol; 8, 7-ketocholesterol. All compounds were derivatized with TMS. Reprinted with permission.

Figure 3. LC-MS performance of the method by S Matysik and co-workers [42]. Selected peaks: **1.** 25-hydroxycholesterol, **2.** 24(S)-hydroxycholesterol, **3.** 27-hydroxycholesterol, **4.** 7β -hydroxycholesterol, **5.** 7α -hydroxycholesterol, **6.** 4β -hydroxycholesterol, **7.** 7-ketocholesterol, **8.** cholestan- 3β , 5α , 6β -triol. All compounds were derivatized with DMG. Reprinted with permission.

Figure 4. The most common derivatizations reaction end products for oxysterol analysis; Girard P, Girard T, Picolinyl ester and N,N-dimethylglycin.

Figure 5. GCxGC-FID performance of the method by Tranchida et al. [51]. The 2D chromatogram is of a commercial sunflower oil («sterol zone»). Compounds are: 1. Cholesterol (methylsterol = DesMe), 2. Brassicasterol (DesMe), 3. Ergosta-5,7,9(11),22-tetraen-3beta-ol (DesMe), 4. Ergosterol (DesMe), 5. 24-methylene-cholesterol (DesMe), 6. Campesterol (DesMe), 7. Campestanol (DesMe), 8. Stigmasterol (DesMe), 9. Ergosta-7-en-3β-ol (DesMe), 10. Clerosterol (DesMe), 11. β-sitosterol (DesMe), 12. Lupeol (dimethylsterol = DiMe), 13. Δ⁵-avenasterol (DesMe), 14. Parkeol (DiMe), 15. β-amyrin (DiMe), 16. Δ⁷-stigmasterol (DesMe), 20. α-amyrin (DiMe), 21. 24-methylene-cycloartanol (DiMe), 22. Citrostadienol (methylsterol = Me), 23. Erythrodiol (Diol). Reprinted with permission.

References

 C.D. DuSell, M. Umetani, P.W. Shaul, D.J. Mangelsdorf, D.P. McDonnell, 27-Hydroxycholesterol Is an Endogenous Selective Estrogen Receptor Modulator, Molecular Endocrinology, 22 (2008) 65-77.
 S. Nachtergaele, D.M. Whalen, L.K. Mydock, Z. Zhao, T. Malinauskas, K. Krishnan, P.W. Ingham, D.F. Covey, C. Siebold, R. Rohatgi, Structure and function of the Smoothened extracellular domain in

vertebrate Hedgehog signaling, eLife, 2 (2013) e01340.

[3] H. Roberg-Larsen, M.F. Strand, S. Krauss, S.R. Wilson, Metabolites in vertebrate Hedgehog signaling, Biochemical and Biophysical Research Communications, 446 (2014) 669-674.

[4] W.J. Griffiths, J. Abdel-Khalik, P.J. Crick, E. Yutuc, Y. Wang, New methods for analysis of oxysterols and related compounds by LC–MS, The Journal of Steroid Biochemistry and Molecular Biology, 162 (2016) 4-26.

[5] R. Narayanaswamy, V. Iyer, P. Khare, M.L. Bodziak, D. Badgett, R. Zivadinov, B. Weinstock-Guttman, T.C. Rideout, M. Ramanathan, R.W. Browne, Simultaneous Determination of Oxysterols, Cholesterol and 25-Hydroxy-Vitamin D3 in Human Plasma by LC-UV-MS, PLoS ONE, 10 (2015) e0123771.

[6] F.D. Porter, D.E. Scherrer, M.H. Lanier, S.J. Langmade, V. Molugu, S.E. Gale, D. Olzeski, R. Sidhu, D.J. Dietzen, R. Fu, C.A. Wassif, N.M. Yanjanin, S.P. Marso, J. House, C. Vite, J.E. Schaffer, D.S. Ory, Cholesterol Oxidation Products Are Sensitive and Specific Blood-Based Biomarkers for Niemann-Pick C1 Disease, Science Translational Medicine, 2 (2010) 56ra81-56ra81.

[7] H. Roberg-Larsen, C. Vesterdal, S.R. Wilson, E. Lundanes, Underivatized oxysterols and nanoLC– ESI-MS: A mismatch, Steroids, 99 (2015) 125-130.

[8] E. Stokvis, H. Rosing, J.H. Beijnen, Stable isotopically labeled internal standards in quantitative bioanalysis using liquid chromatography/mass spectrometry: necessity or not?, Rapid Communications in Mass Spectrometry, 19 (2005) 401-407.

[9] L.L. Jessome, D.A. Volmer, Ion Suppression: A Major Concern in Mass Spectrometry, LCGC North America, 24 (2016) 498-510.

[10] R.K. Beerthuis, J.H. Recourt, Sterol Analysis by Gas Chromatography, Nature, 186 (1960) 372-374.

[11] C.J.W. Brooks, E.C. Horning, J.S. Young, Characterization of sterols by gas chromatography-mass spectrometry of the trimethylsilyl ethers, Lipids, 3 (1968) 391-402.

[12] R.A. Landowne, S.R. Lipsky, The Electron Capture Spectrometry of Haloacetates: A Means of Detecting Ultramicro Quantities of Sterols by Gas Chromatography, Analytical Chemistry, 35 (1963) 532-535.

[13] S. Dzeletovic, O. Breuer, E. Lund, U. Diczfalusy, Determination of Cholesterol Oxidation Products in Human Plasma by Isotope Dilution-Mass Spectrometry, Analytical Biochemistry, 225 (1995) 73-80. [14] F. Kannenberg, J.-R. Nofer, E. Schulte, J. Reunert, T. Marquardt, M. Fobker, Determination of serum cholestane- 3β , 5α , 6β -triol by gas chromatography–mass spectrometry for identification of Niemann-Pick type C (NPC) disease, The Journal of Steroid Biochemistry and Molecular Biology, 169 (2017) 54-60.

[15] V. Cardenia, M.T. Rodriguez-Estrada, E. Baldacci, S. Savioli, G. Lercker, Analysis of cholesterol oxidation products by Fast gas chromatography/mass spectrometry, Journal of Separation Science, 35 (2012) 424-430.

[16] H. Roberg-Larsen, High sensitivity measurements of side-chain oxysterols using liquid chromatography mass spectrometry, Ph D Thesis. Department of Chemistry, University of Oslo, 2015, Table 3, Page 16.

[17] Z. Zhang, D. Li, D.E. Blanchard, S.R. Lear, S.K. Erickson, T.A. Spencer, Key regulatory oxysterols in liver: analysis as Δ4-3-ketone derivatives by HPLC and response to physiological perturbations, Journal of Lipid Research, 42 (2001) 649-658.

[18] Q. Bai, X. Zhang, L. Xu, G. Kakiyama, D. Heuman, A. Sanyal, W.M. Pandak, L. Yin, W. Xie, S. Ren, Oxysterol sulfation by cytosolic sulfotransferase suppresses liver X receptor/sterol regulatory

element binding protein–1c signaling pathway and reduces serum and hepatic lipids in mouse models of nonalcoholic fatty liver disease, Metabolism, 61 (2012) 836-845.

[19] M. Careri, D. Ferretti, P. Manini, M. Musci, Evaluation of particle beam high-performance liquid chromatography–mass spectrometry for analysis of cholesterol oxides, Journal of Chromatography A, 794 (1998) 253-262.

[20] C.F. Poole, N. Lenca, Applications of the solvation parameter model in reversed-phase liquid chromatography, Journal of Chromatography A, 1486 (2017) 2-19.

[21] N.C. van de Merbel, K.J. Bronsema, M.W.J. van Hout, R. Nilsson, H. Sillén, A validated liquid chromatography–tandem mass spectrometry method for the quantitative determination of 4β-hydroxycholesterol in human plasma, Journal of Pharmaceutical and Biomedical Analysis, 55 (2011) 1089-1095.

[22] H. Shan, J. Pang, S. Li, T.B. Chiang, W.K. Wilson, G.J. Schroepfer, Chromatographic behavior of oxygenated derivatives of cholesterol, Steroids, 68 (2003) 221-233.

[23] H. Roberg-Larsen, K. Lund, T. Vehus, N. Solberg, C. Vesterdal, D. Misaghian, P.A. Olsen, S. Krauss, S.R. Wilson, E. Lundanes, Highly automated nano-LC/MS-based approach for thousand cell-scale quantification of side chain-hydroxylated oxysterols, Journal of Lipid Research, 55 (2014) 1531-1536.
[24] H. Roberg-Larsen, K. Lund, K.E. Seterdal, S. Solheim, T. Vehus, N. Solberg, S. Krauss, E. Lundanes, S.R. Wilson, Mass spectrometric detection of 27-hydroxycholesterol in breast cancer exosomes, The Journal of Steroid Biochemistry and Molecular Biology, 169 (2017) 22-28.

[25] H. Roberg-Larsen, M.F. Strand, A. Grimsmo, P.A. Olsen, J.L. Dembinski, F. Rise, E. Lundanes, T. Greibrokk, S. Krauss, S.R. Wilson, High sensitivity measurements of active oxysterols with automated filtration/filter backflush-solid phase extraction-liquid chromatography–mass spectrometry, Journal of Chromatography A, 1255 (2012) 291-297.

[26] I. Burkard, K.M. Rentsch, A. von Eckardstein, Determination of 24S- and 27-hydroxycholesterol in plasma by high-performance liquid chromatography-mass spectrometry, Journal of Lipid Research, 45 (2004) 776-781.

[27] A.E. DeBarber, D. Lütjohann, L. Merkens, R.D. Steiner, Liquid chromatography–tandem mass spectrometry determination of plasma 24S-hydroxycholesterol with chromatographic separation of 25-hydroxycholesterol, Analytical Biochemistry, 381 (2008) 151-153.

[28] E. Lesellier, A simple and interactive column classification for reversed-phase liquid chromatography: The carotenoid test, Part 1: Studied properties, probes, and fundametnals, LCGC North America, 34 (2016) 766-777.

[29] Cartenoid test for classification map of reveres phase columns, Institut de Chimie Organique et Analytique, 23/02/2018, <u>http://www.icoa.fr/classification_c18_columns/index.php</u>

[30] P.J. Crick, L. Beckers, M. Baes, P.P. Van Veldhoven, Y. Wang, W.J. Griffiths, The oxysterol and cholestenoic acid profile of mouse cerebrospinal fluid, Steroids, 99 (2015) 172-177.

[31] P.J. Crick, W.J. Griffiths, J. Zhang, M. Beibel, J. Abdel-Khalik, J. Kuhle, A.W. Sailer, Y. Wang, Reduced Plasma Levels of 25-Hydroxycholesterol and Increased Cerebrospinal Fluid Levels of Bile Acid Precursors in Multiple Sclerosis Patients, Molecular Neurobiology, 54 (2017) 8009-8020.

[32] P.J. Crick, T. William Bentley, J. Abdel-Khalik, I. Matthews, P.T. Clayton, A.A. Morris, B.W. Bigger, C. Zerbinati, L. Tritapepe, L. Iuliano, Y. Wang, W.J. Griffiths, Quantitative Charge-Tags for Sterol and Oxysterol Analysis, Clinical Chemistry, 61 (2015) 400-411.

[33] W.J. Griffiths, J. Abdel-Khalik, P.J. Crick, M. Ogundare, C.H. Shackleton, K. Tuschl, M.K. Kwok, B.W. Bigger, A.A. Morris, A. Honda, L. Xu, N.A. Porter, I. Björkhem, P.T. Clayton, Y. Wang, Sterols and oxysterols in plasma from Smith-Lemli-Opitz syndrome patients, The Journal of Steroid Biochemistry and Molecular Biology, 169 (2017) 77-87.

[34] K. Karu, M. Hornshaw, G. Woffendin, K. Bodin, M. Hamberg, G. Alvelius, J. Sjövall, J. Turton, Y. Wang, W.J. Griffiths, Liquid chromatography-mass spectrometry utilizing multi-stage fragmentation for the identification of oxysterols, Journal of Lipid Research, 48 (2007) 976-987.

[35] A. Honda, T. Miyazaki, T. Ikegami, J. Iwamoto, K. Yamashita, M. Numazawa, Y. Matsuzaki, Highly sensitive and specific analysis of sterol profiles in biological samples by HPLC–ESI–MS/MS, The Journal of Steroid Biochemistry and Molecular Biology, 121 (2010) 556-564.

[36] A. Honda, K. Yamashita, T. Hara, T. Ikegami, T. Miyazaki, M. Shirai, G. Xu, M. Numazawa, Y. Matsuzaki, Highly sensitive quantification of key regulatory oxysterols in biological samples by LC-ESI-MS/MS, Journal of Lipid Research, 50 (2009) 350-357.

[37] T. Ikegami, A. Honda, T. Miyazaki, M. Kohjima, M. Nakamuta, Y. Matsuzaki, Increased serum oxysterol concentrations in patients with chronic hepatitis C virus infection, Biochemical and Biophysical Research Communications, 446 (2014) 736-740.

[38] Y. Xu, Y. Yuan, L. Smith, R. Edom, N. Weng, R. Mamidi, J. Silva, D.C. Evans, H.-K. Lim, LC–ESI-MS/MS quantification of 4β-hydroxycholesterol and cholesterol in plasma samples of limited volume, Journal of Pharmaceutical and Biomedical Analysis, 85 (2013) 145-154.

[39] M.-Q. Huang, W. Lin, W. Wang, W. Zhang, Z. Lin, N. Weng, Quantitation of P450 3A4 endogenous biomarker – 4β -hydroxycholesterol – in human plasma using LC/ESI-MS/MS, Biomedical Chromatography, 28 (2014) 794-801.

[40] R. Hayes, A. Ahmed, T. Edge, H. Zhang, Core–shell particles: Preparation, fundamentals and applications in high performance liquid chromatography, Journal of Chromatography A, 1357 (2014) 36-52.

[41] J.G. McDonald, D.D. Smith, A.R. Stiles, D.W. Russell, A comprehensive method for extraction and quantitative analysis of sterols and secosteroids from human plasma, Journal of Lipid Research, 53 (2012) 1399-1409.

[42] Z. Pataj, G. Liebisch, G. Schmitz, S. Matysik, Quantification of oxysterols in human plasma and red blood cells by liquid chromatography high-resolution tandem mass spectrometry, Journal of Chromatography A, 1439 (2016) 82-88.

[43] D.S. Bell, D. Shollenberger, H. Cramer, Evaluation of retention and selectivity using biphenyl stationary phases, LCGC North America, 35 (2017) 360-365.

[44] X. Jiang, R. Sidhu, F.D. Porter, N.M. Yanjanin, A.O. Speak, D.T. te Vruchte, F.M. Platt, H. Fujiwara, D.E. Scherrer, J. Zhang, D.J. Dietzen, J.E. Schaffer, D.S. Ory, A sensitive and specific LC-MS/MS method for rapid diagnosis of Niemann-Pick C1 disease from human plasma, Journal of Lipid Research, 52 (2011) 1435-1445.

[45] S. Boenzi, F. Deodato, R. Taurisano, D. Martinelli, D. Verrigni, R. Carrozzo, E. Bertini, A. Pastore, C. Dionisi-Vici, D.W. Johnson, A new simple and rapid LC–ESI-MS/MS method for quantification of plasma oxysterols as dimethylaminobutyrate esters. Its successful use for the diagnosis of Niemann–Pick type C disease, Clinica Chimica Acta, 437 (2014) 93-100.

[46] M. Romanello, S. Zampieri, N. Bortolotti, L. Deroma, A. Sechi, A. Fiumara, R. Parini, B. Borroni, F. Brancati, A. Bruni, C.V. Russo, A. Bordugo, B. Bembi, A. Dardis, Comprehensive Evaluation of Plasma 7-Ketocholesterol and Cholestan- 3β , 5α , 6β -Triol in an Italian Cohort of Patients Affected by Niemann-Pick Disease due to NPC1 and SMPD1 Mutations, Clinica Chimica Acta, 455 (2016) 39-45.

[47] S. Adhikari, Y. Xia, Thiyl Radical-Based Charge Tagging Enables Sterol Quantitation via Mass Spectrometry, Analytical Chemistry, 89 (2017) 12631-12635.

[48] E. Cha, K.M. Lee, K.D. Park, K.S. Park, K.-W. Lee, S.-M. Kim, J. Lee, Hydroxycholesterol Levels in the Serum and Cerebrospinal Fluid of Patients with Neuromyelitis Optica Revealed by LC-Ag+CIS/MS/MS and LC-ESI/MS/MS with Picolinic Derivatization: Increased Levels and Association with Disability during Acute Attack, PLOS ONE, 11 (2016) e0167819.

[49] T. Vehus, H. Roberg-Larsen, J. Waaler, S. Aslaksen, S. Krauss, S.R. Wilson, E. Lundanes, Versatile, sensitive liquid chromatography mass spectrometry – Implementation of 10 μm OT columns suitable for small molecules, peptides and proteins, Scientific Reports, 6 (2016) 37507.

[50] P.Q. Tranchida, F.A. Franchina, P. Dugo, L. Mondello, Comprehensive two-dimensional gas chromatography-mass spectrometry: Recent evolution and current trends, Mass Spectrometry Reviews, 35 (2016) 524-534.

[51] P.Q. Tranchida, S. Salivo, F.A. Franchina, I. Bonaccorsi, P. Dugo, L. Mondello, Qualitative and quantitative analysis of the unsaponifiable fraction of vegetable oils by using comprehensive 2D GC with dual MS/FID detection, Analytical and Bioanalytical Chemistry, 405 (2013) 4655-4663.

[52] M. Holčapek, M. Ovčačíková, M. Lísa, E. Cífková, T. Hájek, Continuous comprehensive twodimensional liquid chromatography–electrospray ionization mass spectrometry of complex lipidomic samples, Analytical and Bioanalytical Chemistry, 407 (2015) 5033-5043.

[53] A. Periat, J. Boccard, J.-L. Veuthey, S. Rudaz, D. Guillarme, Systematic comparison of sensitivity between hydrophilic interaction liquid chromatography and reversed phase liquid chromatography coupled with mass spectrometry, Journal of Chromatography A, 1312 (2013) 49-57.

[54] M. Gregus, H. Roberg-Larsen, E. Lundanes, F. Foret, P. Kuban, S.R. Wilson, Non-aqueous capillary electrophoretic separation of cholesterol and 25-hydroxycholesterol after derivatization with Girard P reagent, Chemistry and Physics of Lipids, 207 (2017) 87-91.

[55] T. Hara, S. Futagami, S. Eeltink, W. De Malsche, G.V. Baron, G. Desmet, Very high efficiency porous silica layer open-tubular capillary columns produced via in-column sol–gel processing, Analytical chemistry, 88 (2016) 10158-10166.

[56] Y.M. Choo, M.H. Ng, A.N. Ma, C.H. Chuah, M.A. Hashim, Application of supercritical fluid chromatography in the quantitative analysis of minor components (carotenes, vitamin E, sterols, and squalene) from palm oil, Lipids, 40 (2005) 429-432.

[57] K.-H. Storbeck, L. Gilligan, C. Jenkinson, E.S. Baranowski, J.L. Quanson, W. Arlt, A.E. Taylor, The utility of ultra-high performance supercritical fluid chromatography–tandem mass spectrometry (UHPSFC-MS/MS) for clinically relevant steroid analysis, Journal of Chromatography B, 1085 (2018) 36-41.

[58] S. Ayciriex, A. Regazzetti, M. Gaudin, E. Prost, D. Dargère, F. Massicot, N. Auzeil, O. Laprévote, Development of a novel method for quantification of sterols and oxysterols by UPLC-ESI-HRMS: application to a neuroinflammation rat model, Analytical and Bioanalytical Chemistry, 404 (2012) 3049-3059.

[59] G. Klinke, M. Rohrbach, R. Giugliani, P. Burda, M.R. Baumgartner, C. Tran, M. Gautschi, D. Mathis, M. Hersberger, LC-MS/MS based assay and reference intervals in children and adolescents for oxysterols elevated in Niemann–Pick diseases, Clinical Biochemistry, 48 (2015) 596-602.

[60] A. Honda, K. Yamashita, H. Miyazaki, M. Shirai, T. Ikegami, G. Xu, M. Numazawa, T. Hara, Y. Matsuzaki, Highly sensitive analysis of sterol profiles in human serum by LC-ESI-MS/MS, Journal of Lipid Research, 49 (2008) 2063-2073.

[61] A.K. Goodenough, J.M. Onorato, Z. Ouyang, S. Chang, A.D. Rodrigues, S. Kasichayanula, S.-P. Huang, W. Turley, R. Burrell, M. Bifano, M. Jemal, F. LaCreta, A. Tymiak, D. Wang-Iverson, Quantification of 4-Beta-Hydroxycholesterol in Human Plasma Using Automated Sample Preparation and LC-ESI-MS/MS Analysis, Chemical Research in Toxicology, 24 (2011) 1575-1585.

[62] P.J. Crick, J. Aponte, T.W. Bentley, I. Matthews, Y. Wang, W.J. Griffiths, Evaluation of novel derivatisation reagents for the analysis of oxysterols, Biochemical and Biophysical Research Communications, 446 (2014) 756-761.

[63] L. Baila-Rueda, A. Cenarro, M. Cofan, I. Orera, S. Barcelo-Batllori, M. Pocovi, E. Ros, F. Civeira, C. Nerin, C. Domeno, Simultaneous determination of oxysterols, phytosterols and cholesterol precursors by high performance liquid chromatography tandem mass spectrometry in human serum, Analytical Methods, 5 (2013) 2249-2257.

[64] A. Gorassini, G. Verardo, S.-C. Fregolent, R. Bortolomeazzi, Rapid determination of cholesterol oxidation products in milk powder based products by reversed phase SPE and HPLC-APCI-MS/MS, Food Chemistry, 230 (2017) 604-610.

[65] M. Hasan, W. Siegmund, S. Oswald, Rapid LC–MS/MS method for the determination of 4hydroxycholesterol/cholesterol ratio in serum as endogenous biomarker for CYP3A activity in human and foals, Journal of Chromatography B, 1033-1034 (2016) 193-199.

[66] C. Helmschrodt, S. Becker, J. Schröter, M. Hecht, G. Aust, J. Thiery, U. Ceglarek, Fast LC–MS/MS analysis of free oxysterols derived from reactive oxygen species in human plasma and carotid plaque, Clinica Chimica Acta, 425 (2013) 3-8.

[67] R. Karuna, A. von Eckardstein, K.M. Rentsch, Dopant assisted-atmospheric pressure photoionization (DA-APPI) liquid chromatography–mass spectrometry for the quantification of 27-hydroxycholesterol in plasma, Journal of Chromatography B, 877 (2009) 261-268.

[68] P. Manini, Andreoli, R., Careri, M., Elviri, L. and Musci, M., Atmospheric pressure chemical ionization liquid chromatography/mass spectrometry in cholesterol oxide determination and characterization., Rapid Commun. Mass Spectrom., 12 (1998) 883-889.

[69] I. Mendiara, C. Domeño, C. Nerín, A.M. Geurts, J. Osada, R. Martínez-Beamonte, Determination of total plasma oxysterols by enzymatic hydrolysis, solid phase extraction and liquid chromatography coupled to mass-spectrometry, Journal of pharmaceutical and biomedical analysis, 150 (2018) 396-405.

[70] K. Raith, C. Brenner, H. Farwanah, G. Müller, K. Eder, R.H.H. Neubert, A new LC/APCI-MS method for the determination of cholesterol oxidation products in food, Journal of Chromatography A, 1067 (2005) 207-211.

[71] E. Razzazi-Fazeli, S. Kleineisen, W. Luf, Determination of cholesterol oxides in processed food using high-performance liquid chromatography–mass spectrometry with atmospheric pressure chemical ionisation, Journal of Chromatography A, 896 (2000) 321-334.

[72] T. Saldanha, A.C.H.F. Sawaya, M.N. Eberlin, N. Bragagnolo, HPLC Separation and Determination of 12 Cholesterol Oxidation Products in Fish: Comparative Study of RI, UV, and APCI-MS Detectors, Journal of Agricultural and Food Chemistry, 54 (2006) 4107-4113.

[73] I.H.K. Dias, I. Milic, G.Y.H. Lip, A. Devitt, M.C. Polidori, H.R. Griffiths, Simvastatin reduces circulating oxysterol levels in men with hypercholesterolemia, Redox Biology, (2018).





Highlights

- ➤ Gas chromatography for oxysterols is well established
- > More diversity in liquid chromatography approaches
- > Approaches to enhancing speed and selectivity with LC are presented
- > Few differences in chromatography between native and derivatized oxysterols
- > Alternative approaches to separation of oxysterols exist, but are little explored